

ONEPro

MOBILE MESH NEBULIZER



Always with you

Very easy to carry and use in any travel condition thanks to its small size.



Full Green Energy Saving

No internal battery to recharge. The device is powered by the USB devices you use daily.



Quiet

The peculiar design and the small size allow to apply it to the face and to use it also hands-free.



USB Powered

The device can be powered through any USB device.

INCLUDED ACCESSORIES

ADULT MASK + CHILD MASK

with elastic strings

MOUTHPIECE + NOSEPIECE

non-invasive

POWER SUPPLY

2 USB ADAPTER

USB-C and USB-Lightning

PRODUCT CODE RE 300912/01

KIWI+

MOBILE MESH NEBULIZER



Fast & efficient

Second generation mesh nebulizer, 8ml capacity, takes only 6 minutes with 2,5 ml of drug (0,30ml/min).



Long lasting batteries

The lithium batteries last up to 1100 hours of operation or 300 charging cycles, making it the ideal device for travelling.



Always with you

Tiny and easily transportable thanks to its small size. It can be operated by lithium battery and easily stored.



Baby friendly

Easy to hold thanks to its small size and pleasant to hear thanks to its low noise level (38,2dB(A)),



High-quality accessories

Equipped with a wide range of accessories with comfortable soft touch finishings.



Autocleaning function

Reduces the risk of annoying clogging of the mesh caused by improper cleaning at the end of each use.



Excellent nebulization

Fine particle size of 1,51 micron, 93% of which is <5 micron and <0,1ml residual volume.

INCLUDED ACCESSORIES

VIBRATING MESH CHAMBER

2nd generation

ADULT MASK + CHILD MASK

with soft touch finishing

MOUTHPIECE + NOSEPIECE

non-invasive

CARRYING CASE

for safe storage and transport

AC/DC UNIVERSAL POWER ADAPTER

100-240V 50-60Hz

LITHIUM BATTERY

PRODUCT CODE RE 300911